

Lipo laser Vans3




1. Application

- a. intensive physicaallipolysis to remove fat;
- b. excess fat cell melted;
- c. body slimming, cellulite reduction;
- d. smooth fatigue;
- e. remove obstruction from channels and collaterals;
- f. promote and accelerate the body's metabolism;

2. Specification

Model	Vans3
Laser wavelength	650nm
Laser type	Diode
Numbers of treatment probe	10 (Eight big ones and two small ones)
Numbers of diode laser	34
Output Power	40w
Timer	30 minutes max
Output Frequency	1 hz to 1000hz adjustable
Cooling Requirements	Air cooling
Working temperature	10-40 degree
Electrical Source	10V/120V-220V/240V

3. menu



4. real photo

